LIFEWAVE®

SP6®
Complete

Instructions

Anwendung

Instrucciones

Mode d'emploi

Leiðbeiningar

Istruzioni

Instruções

Instructions for Use


SP6 Complete should be used 5 days per week, wearing a new patch each day. Suggested use for week 1 shown in diagrams. For maximum results during week 2 and thereafter, rotate the patch only between the points on your left leg that give you the most control of your appetite.

Apply the patches to clean, dry skin in the morning. Patches may be worn for up to 12 hours before discarding. Keep well hydrated while using this product. Do not use more than one patch at a time on the body.

STEP 1

First, locate the acupressure point you wish to stimulate.

STEP 2

Next, remove a plastic bead from the bag. Remove the backing from the adhesive plastic patch included. Place the plastic bead in the center of the patch on the adhesive side of the patch.

STEP 3

Lastly, apply the bead and patch to the acupressure point, smoothing the patch on so that it is properly secured to the skin. There should be a mild pressure on the acupressure point.

Warnings: Remove immediately if you feel discomfort or skin irritation occurs. Do not reuse patch once removed from the skin. For external use only. Do not ingest. Do not use on wounds or damaged skin. Ask a health professional before using if you have a health condition, any questions or concerns about your health. Do not use if pregnant or nursing. This booklet is no way intended to be medical advice. Individual results will vary.

DAY 1 - SPLEEN 6

Located about four finger widths above the inside ankle bone of the LEFT leg, on the tibia. This point has given the best results.

DAY 2 - STOMACH 36

Located in the depression 4 finger widths below the LEFT knee cap and 1 finger widths outside the bone. This point gives good results.

DAY 3 - KIDNEY 3

Located at the depression on the inside of the LEFT ankle, just behind the ankle bone

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Manufactured by: LifeWave Europe, Ltd. | Raheen Ind. Est., Athenry, Co. Galway, IRL
European Authorized Representative: Obelis S.A. | Bd. Général Wahis, 53, 1030 Brussels,
Hergestellt von: LifeWave Europe, Ltd. | Raheen Ind. Est., Athenry, Co. Galway, IRL
Autorisierte europäische Vertretung: Obelis S.A. | Bd. Général Wahis, 53, 1030 Brüssel,
Fabricado por: LifeWave Europe, Ltd. | Raheen Ind. Est., Athenry, Co. Galway, IRL
Representante autorizado en Europa: Obelis S.A. | Bd. Général Wahis, 53, 1030 Bruselas,