

Using JDICOM

Introduction

JDICOM is a summary of DICOM applications which are very useful for troubleshooting. Our best appreciation to Mr. Gunter Zeilinger which created this great tool.

If you have any questions or problems, do not hesitate to contact us. The e-mail addresses are

ccc@med.siemens.de or klaus.gartner@med.siemens.de.

The description is valid for the actual version of JDICOM 1.7.5. It was installed on the following operating systems Win98, WinNT, Win2000, Linux (RedHat 6.1 and Suse 7.0). Required is a java2 runtime environment for the used hardware architecture.

Functionality of contained applications

ImageServer

Supports Image related DICOM Services Classes:

Receives images over network and stores them on (local) disk

(= SCP of Storage SOP classes + FSC+FSU of Media Storage SOP classes).

Dispatches received query/retrieve requests by transmitting information about stored images or images themselves to specified destination.

(= SCP of Query/Retrieve SOP classes + FSR of Media Storage SOP classes).

Notifies specified peer application about received images

(= SCU of Study Content Notification SOP class).

Sends requested commitments back about received images.

(=SCU of Storage Commitment SOP class)

StorageSCU

Load images from (local) disk and send them over network.

(= SCU of Storage SOP classes).

Receives commitments back about transmitted images.

(=SCP of Storage Commitment SOP class)

QueryRetrieveSCU

Queries image archives and controls remote retrieve of images to specified destination.

(= SCU of Query/Retrieve SOP classes).

RisServer

Supports RIS related DICOM Services Classes:

Receives worklist requests over the network

(= SCP of Basic Worklist SOP class).

Receives notifications about the (image) content of studies over network

(= SCP of Study Content Notification SOP class).

Receives performed procedure steps about the patient from the network

(= SCU of Performed Procedure Step SOP class).

ModalitySCU

Supports RIS related DICOM Services Classes:

Sends worklist requests over the network

(= SCU of Basic Worklist SOP class).

Receives notifications about the (image) content of studies over network

(= SCP of Study Content Notification SOP class).

Sends performed procedure steps about the patient over the network to the RisServer.

(= SCP of Performed Procedure Step SOP class).

PrintServer

Supports Basic Print related DICOM Services Classes:
Receives printable images over network
(= SCP Basic Print SOP class).

PrintSCU (not yet described)

EditDicomObject

Administration tool for DICOM File to:
Display and modify the DICOM files tags.

EditDicomDir

Administration tool for DICOM File-sets to:
Display and modify the DICOMDIR file containing the Media Storage Directory information (= FSC+FSU of Media Storage SOP classes). Delete images from the File-sets, by deleting the image files and updating the DICOMDIR file accordingly. Import images into the File-sets, by updating the DICOMDIR file with Directory information extracted from the new added images.

Validate

Checks DICOM files - containing Composite Information Objects - for its DICOM conformance.
Lists attribute names and values of contained Data Elements.

Getting started

Installation

The JDicom applications runs on Java Platform 2. So you can either download the JDicom version which already includes the Java virtual machine from the following location :

<http://www.tiani.com/JDicom/install/>

Use the link for Windows "[Windows including a JVM](#)" and for Linux "[Linux including a JVM](#)" in order to download the corresponding JDicom package.

If the file is successfully downloaded use following the installation instruction for Windows:



To install this application on a Windows system:

1. Save [jinstall.exe](#) in a temporary directory on your system as a **binary** file. Important: To download a binary file and save it under the correct name, you usually have to right click on the filename's jump point and select Save from the popup menu.
2. Verify that the filename you downloaded is **jinstall.exe**.
If your browser saved it with any other name, then rename it.
3. Double-click on the downloaded copy of **jinstall.exe** to install the new application.

for Linux:




To install this application on a Windows system:

1. Save [jinstall.exe](#) in a temporary directory on your system as a **binary** file. Important: To download a binary file and save it under the correct name, you usually have to right click on the filename's jump point and select Save from the popup menu.
2. Verify that the filename you downloaded is **jinstall.exe**.
If your browser saved it with any other name, then rename it.
3. Double-click on the downloaded copy of **jinstall.exe** to install the new application.

to install JDicom on your PC.

After the successful installation of JDicom you find this menu structure in "START -> PROGRAMS"



The important entries are:

EditDicomDir
EditDicomObjects
ImageServer
ImageServerGUI
ModalitySCU
PrintSCU
PrintServer
PrintServerGUI
QueryRetrieveSCU
RisServer
RisServerGUI
StorageSCU
Validate

The files are installed by using the default settings in the directory: "C:\Program Files\jdicom"
Automatically subfolders are created in the default installation directory. In the subdirectory applet you can find the available applet version of JDicom. These applets are running in a Web Browser which supports Java.


Quick tour

The following DICOM services are not yet described in this document:


- ◆ Storage Commitment SCU/SCP
- ◆ Query/Retrieve SCP
- ◆ Print SCU


StorageSCU

With the JDICOM application StorageSCU you can simulate the behaviour of the DICOM service Storage as Service Class User. An Association can be opened, a DICOM Verification can be made and DICOM images can be transferred to a remote DICOM destination.



The DICOM application StorageSCU runs with a graphical user interface (GUI). Select the entry

 **StorageSCU** in the JDicom menu.



Here are the minimum entries which could be modified before you can use StorageSCU:

Host :

Enter here the IP address or hostname of the StorageSCP (DICOM partner for receiving the image)

If a host name is used, an entry in the hosts file is required or DNS has to be enabled.

Port :

Portnumber of the StorageSCP (DICOM partner for receiving the image)

Called Title :

DICOM AET of the StorageSCP (DICOM partner for receiving the image)

Calling Title :

Own DICOM AET

Transfer Syntax :

This are the possible values:

ImplicitVRLittleEndian

ExplicitVRLittleEndian

ExplicitVRBigEndian

JPEGLossless

JPEGBaseline (lossy for 8 bit images)

If you are using a wrong entry for the TransferSyntax a message is displayed in the log tab when you try to open the association with "**Connect**":

TransferSyntax: java.lang.IllegalArgumentException: not one of {ImplicitVRLittleEndian,ExplicitVRLittleEndian,ExplicitVRBigEndian,JPEGLossless,JPEGBaseline}

With the following properties it is possible to change the Patient Name and Patient ID during the image transfer to the new values:

Anonymize: **false** (no change is done)
true (change the name and ID)

Anonymize.NewName: new PatientName

Anonymize.NewID : new PatientID

Verbose :

Specifies the log level. The value which can be used starts with

"0" no output in the log tab window

"1"

"2" decode of DIMSE (DICOM Message Service Element)

"3" decode of data elements

Note: If you change an entry in the Props tab you have to confirm it by pressing the <ENTER> key on your keyboard. If you have changed a property value you have to CLOSE/OPEN the association in order to activate the made changes.

In this example the DICOM partner for receiving images (StorageSCP) has the parameters:

Hostname: **146.254.123.173**

Portnumber of the StoreSCP process: **104**

AET of the StoreSCP process: **AN_TDSIM_1**

The local (own) DICOM AET of the StorageSCU JDICOM program: **Storage_SCU**

Transfer Syntax: **ImplicitVRLittleEndian**

Anonymize: **false**

Verbose mode: **3**

Start the Association by clicking the "**Connect**" button.

With "**Echo**" you start a Verification (C_Echo).

With "**Send**" only one selected image can be transferred.


With "**Send***" multiple images can be transferred.

With "**Cancel**" the started image transfer can be interrupted.

Close the Association by clicking the "**Release**" button.

ImageServer


With the JDICOM application ImageServer you can simulate the behaviour of the DICOM service Storage as Service Class Provider. DICOM Verification is supported, DICOM images can be received and stored into a specified directory. The ImageServer simulates also the DICOM services Query/Retrieve and Storage Commitment as Provider.



The DICOM application ImageServer can be run either with a graphical user interface (GUI) or in command line mode.

To run it with a GUI select  **ImageServerGUI** in the JDicom menu.

To run it in command line mode select  **ImageServer** from the JDicom menu (recommended only for advanced users).



Fileset.TransferSyntax:
Used Transfer Syntax for storing the images to harddisc. Possible values:

- AcrNema**
- ImplicitVRLittleEndian**
- ExplicitVRLittleEndian**
- ExplicitVRBigEndian**
- JPEG_14_SelectionValue1** (lossless)

In the ImageServerGUI you find three tab cards.
Props: tabcard for changing the properties of the Imageserver
AET's: tabcard only required for using the Retrieve and Storage Commitment function (destinations entries for it: AET, IP address and Port)
Log: tabcard for the log output area

The following entries in the Props tabcard are the minimum once which should be modified before you can use the ImageServer:

Port:
Own Portnumber of the ImageServer (for DICOM services StoreSCP and Q/R SCP)

Called Title:
Own DICOM AET of the ImageServer (for DICOM services StoreSCP and Q/R SCP)

Calling Title:
If no entry was made here, all AET's from the StoreSCU's and Q/R SCU's are accepted. If you want to restrict the access to the Imageserver you can put in here the allowed AET's. The AET's should be separated by comma.

Verification:
Specifies whether the Imageserver accepts a Verification DIMSE (C_Echo). Possible values are:
true Verification is supported
false Verification is not supported

Fileset:
This value should be set to "**true**", in order to save the received images in the fileset.path directory. Possible values are true and false.

Fileset.path:
Existing directory for storing the received images. If this directory doesn't exist on your PC you will get an error message after starting your ImageServer.
Example for the entry: **d:/dicom/JDicom/image**
You can also use the \ to separate the subdirectories.

Verbose:
Specifies the log level. The value which can be used starts with
"0" no output in the log tab window
"1"
"2" decode of DIMSE (DICOM Message Service Element)
"3" decode of data elements

If you are using an wrong entry for the Fileset.TransferSyntax a message is displayed in the log tab when you try to start the ImageServer with "**Start**":

```
java.lang.IllegalArgumentException: Fileset.TransferSyntax:
java.lang.IllegalArgumentException: not one of
{AcrNema, ImplicitVRLittleEndian, ExplicitVRLittleEndian, ExplicitVRBigEndian,
JPEG_14_SelectionValu1}
```

In this example the following properties are used for the ImageServer:

	Property	Value
Portnumber of the ImageServer process:	Port	104
AET of the ImageServer process:	CalledTitles	Storage_SCP
AET's accepted StoreSCU's:	CallingTitles	
Save images on hard disc enabled:	Fileset	true
Path on the hard disc:	Fileset.Path	d:/dicom/JDicom/image
Image store syntax	Fileset.TransferSyntax:	ImplicitVRLittleEndian
Mode of the output log:	Verbose	1

Note: If you change an entry in the Props tab you have to confirm it by pressing the <ENTER> key on your keyboard. If you have changed a property value you have to STOP/START the JDICOM ImageServer in order to activate the made changes.


Save the made changes in the property tab card by pressing the "**Apply**" button. The property values are than saved in the properties file for the ImageServer (ImageServer.properties).

Start the ImageServer by clicking the "**Start**" button.

With "**Stop**" the ImageServer will be stopped.

RisServer


With the JDICOM application RisServer you can simulate the behaviour of the RIS or HIS system. The JDICOM RisServer supports the DICOM Services Modality Worklist, Modality Performed Procedure Step and Study Content Notification. The JDICOM RisServer is for all mentioned services Service Class Provider. The JDICOM RisServer supports DICOM Verification and can be used to generate worklist entries.



The DICOM application RisServer can be run either with a graphical user interface (GUI) or in command line mode.

To run it with a GUI select  **RisServerGUI** in the JDicom menu.

To run it in command line mode select  **RisServer** from the JDicom menu (recommended only for advanced users).



In the RisServerGUI you find three tab cards.
Props : tabcard for for changing the properties of the RisServer
Result : tabcard for creating/listing worklist entries
Log : tabcard for the log output area

The following entries in the Props tabcard are the minimum once which should be modified before you can use the ImageServer:

Port :
 Own Portnumber of the RisServer (for DICOM Basic Worklist SCP)

Called Title :
 Own DICOM AET of the RisServer (for DICOM Basic Worklist SCP)

Calling Title :
 If no entry was made here, all AET's from the Basic Worklist SCU's are accepted. If you want to restrict the access to the Imageserver you can put in here the allowed AET's. The AET's should be separated by comma.

Verification :
 Specifies wether the Rissserver accepts a Verification DIMSE (C_Echo). Possible values are:
true Verification is supported
false Verification is not supported

Repository.path :
 Existing directory for storing the created WL entries. If this directory doesn't exist on your PC you will get an error message after starting your RisServer. Example for the entry : **d:/dicom/jdicom/rissserver**
 You can also use the \ to separate the subdirectories.

Verbose :
 Specifies the log level. The value which can be used starts with
"0" no output in the log tab window
"1"
"2" decode of DIMSE (DICOM Message Service Element)
"3" decode of data elements

ModalityWorklist :
 Specifies whether the Rissserver acts as a Basic Worklist SCP. Possible values are:
true Basic Worklist support is enabled.
false Basic Worklist support is disenabled

ModalityPPS :
 Specifies whether the Rissserver acts as a Modality Performed Procedure Step SCP. Possible values are:
true ModalityPPS support is enabled.
false ModalityPPS support is disenabled

StudyContentNotification :
 Specifies whether the Rissserver acts as a Study Content Notification SCP. Possible values are:
true Study Content Notification support is enabled.
false Study Content Notification support is disenabled

Save the made changes in the property tab card by pressing the **"Save"** button. The property values are than saved in the properties file for the RisServer (RisServer.properties). You have to do it once after the installation because this file doesn't exist directly after the installation. You can choose a different name for the RisServer.properties file. With **"load"** you can load a different named RisServer.properties file which you have created before.

In this example the following properties are used for the RisServer:

	Property	Value
Portnumber of the RisServer process:	Port	4460
AET of the RisServer process:	CalledTitles	Worklist_SCP
AET's accepted Worklist SCU's:	CallingTitles	
Path on the hard disc:	Repository.Path	d:/dicom/jdicom/worklist
Mode of the output log:	Verbose	1


Note: If you change an entry in the Props tab you have to confirm it by pressing the <ENTER> key on your keyboard. If you have changed a property value you have to STOP/START the JDICOM RisServer in order to activate the made changes.

Start the RisServer by clicking the "**Start**" button.

With "**Stop**" the RisServer will be stopped.

Creation of a new Worklist entry:

- ◆ Select the "**Result**" tab card in the RisServerGUI



Double click the folder icon Repository and the tree with the following entries is expanded (see image on the right side):

ModalityWorklist:



Select entry for generation of Basic Worklist data sets.


ModalityPPS:

Select entry for displaying received Modality Performed Procedure Steps data sets.

StudyContentNotification:

Select entry for generation of Study Content Notification data sets.

In order to generate an Basic Worklist data set entry select "ModalityWorklist" in the Repository tree and press the  button. With the  button an data set entry can be deleted.

The  button you can refresh the Repository tree.


Automatically a worklist entry is generated and saved. It has default values for each DICOM tag which are generated.

The values can be modified under the RisServerGUI platform. For adding additional WL tag's, select the appropriate Worklist entry in the Repository/ModalityWorklist tree

and click on the  button.

Modify the values of the corresponding tag's and save the made modifications by clicking on

the  button.



Important worklist tags which should match for the worklist request from the modality (WL SCU):

Modality:

Insert here your modality abbreviation (e.g. MR, CT, AX, US, OT, ...) [in this example: **MR**]

ScheduledStationAETitle:

Insert here the AET for which the worklist entry is. Only for this AET it is possible to get this specific worklist entry. Is this tag empty all Worklist SCU's from the same modality can retrieve this Worklist entry. [in this example: **HRI_TDSIM_1**]

ScheduledProcedureStepStartDate:

Default date when the WL entry is generated. Change it when the WL entry is not generated today to the date of today.

Note: If you change a tag value you have to confirm it by pressing the <ENTER> key on your keyboard.

To check if your made changes are OK, select the folder ModalityWorklist and then again the WL entry. Now the modified WL file should be loaded into the RisServerGUI.


Start the RisServer by clicking the "**Start**" button if not already done. Check in the log tab card whether the RisServer is successfully started or not.

The messages "Waiting for invocations from clients..." indicates an successful start of the RisServer.

If you get the following message "java.io.IOException: d:\dicom\jdicom\risserv224 is not a writeable directory!" the RisServer is not started. Check the entry for the **Repository.Path** whether this directory exists on your hard disc.

ModalitySCU (WL SCU)


With the JDICOM application ModalitySCU you can simulate the behaviour of the DICOM services Modality Worklist and Modality Performed Procedure Step, both as Service Class User. An Association can be opened, a DICOM Verification can be made and Worklist information can be retrieved. From the retrieved Worklist a performed procedure step can be created and transferred to the RisSystem. In this case the RisSystem must support MPPS.



The DICOM application ModalitySCU runs with a graphical user interface (GUI). Select the entry



in the JDicom menu.



Here are the minimum entries which should be modified in the Property tabcard before you can use ModalitySCU:

Host :

Enter here the IP address or hostname of the RisServer (Worklist SCP)

Port :

Portnumber of the Worklist SCP

Called Title :

DICOM AET of the Worklist SCP (RisServer or HisServer)

Calling Title :

Own DICOM AET

Verbose :

Specifies the log level. The value which can be used starts with

"0" no output in the log tab window

"1"

"2" decode of DIMSE (DICOM Message Service Element)

"3" decode of data elements

There are three more tabcards available:

Filter tabcard for setting the requested WL tags

Result tabcard shows the received worklist entries

Log tabcard shows the log of the DICOM session

Note: If you change an entry in the Props tab you have to confirm it by pressing the <ENTER> key on your keyboard. If you have changed a property value you have to CLOSE/OPEN the association in order to activate the made changes.

In this example the DICOM RisServer (Worklist SCP) has the parameters:

Hostname: **146.254.141.221**

Portnumber of the Worklist SCP process: **4460**

AET of the Worklist SCP process: **Worklist_SCP**

The local (own) DICOM AET of the ModalitySCU JDICOM program: **Worklist_SCU**

Verbose mode: **1**


Start the Association by clicking the "**Open**" button.

With "**Echo**" you start a Verification (C_Echo).

With "**Query WL**" request the worklist entries from the RisServer (C-FIND).

With "**Close**" you can close the association to the RisServer.

- Setting filter in the JDICOM ModalitySCU:




For setting a WL query filter select the Filter tabcard and choose the tag group by selectin one of the tabcards:

- SPS (Scheduled Procedure Step)
- PROCEDURE
- SERVICE REQUEST
- PATIENT
- VISIT

The marked (with ✓) tags are default set. If you want to add tags just mark them. This is done by clicking on . Removing the tag from the WL query just click on the mark ✓ and the is empty.


- Displaying the WL Query result in JDICOM ModalitySCU:


In the Worklist tree you can find the received worklist data sets



Function of the  buttons are:



with this button you can generate an PPS (Performed Procedure Step) entry for the selected worklist entry. In order to do it select first the desired selected worklist entry and click  button. The PPS data set is stored under the folder "*Performed Procedure Steps*".



with this button you can remove an selected data set entry from the folder "*Worklist*" or "*Performed Procedure Steps*".




with this button you can add PPS DICOM tags from the selected PPS entry in the PPS folder.



with this button you can save the made changes in the selected PPS entry.

PrintServerDummy (PrintSCP)


With the JDICOM application PrintServerDummy you can simulate the behaviour of a DICOM Camera. DICOM Verification is supported and a hardcopy file can be received. This file is stored into a specified directory and can be displayed by any DICOM viewer.



The JDICOM application PrintServerDummy can be run either with a graphical user interface (GUI) or in command line mode.

To run it with a GUI select  **PrintServerGUI** in the JDicom menu.

To run it in command line mode select  **PrintServer** from the JDicom menu (recommended only for advanced users).



In the PrintServerDummy GUI you find two tab cards.

Props : tabcard for for changing the properties of the PrintServerDummy.

Log : tabcard for the log output area

The following entries in the Props tabcard are the minimum once which should be modified before you can use the PrintServerDummy.

Port :

Own Portnumber of the PrintServerDummy (for DICOM service Basic Print SCP)

Called Title :

Own DICOM AET of the PrintServerDummy (for DICOM servicesBasic Print SCP)

Calling Title :

If no entry was made here, all AET's from the PrintSCU's are accepted. If you want to restrict the access to the PrintServerDummy you can put in here the allowed AET's. The AET's should be separated by comma.

Verbose :

Specifies the log level. The value which can be used starts with

"0" no output in the log tab window

"1"

"2" decode of DIMSE (DICOM Message Service Element)

"3" decode of data elements

In this example the following properties are used for the PrintServerDummy:

	Property	Value
Portnumber of the PrintServerDummy process:	Port	4450
AET of the PrintServerDummy process:	CalledTitles	P_PC_Klaus
AET's accepted PrintSCU's:	CallingTitles	
Mode of the output log:	Verbose	1

Note: If you change an entry in the Props tab you have to confirm it by pressing the <ENTER> key on your keyboard. If you have changed a property value you have to STOP/START the JDICOM PrintServerDummy in order to activate the made changes.

Save the made changes in the property tab card by pressing the "**Apply**" button. The property values are than saved in the properties file for the PrintServerDummy (PrintServerDummy.properties).

Start the PrintServerDummy by clicking the "**Start**" button.

With "**Stop**" the PrintServerDummy will be stopped.

It is possible to save the received DICOM Basic Print data to hard disc. Therefore you have to modify the following entries in the property tab card.

Property	Value
StoreHC	true
StoreHC.Path	Existing path on PC

In this example the values for this store functionality are :

	Property	Value
Enable hardcopy store to hard disc:	StoreHC	true
Directory for the stored hardcopy files:	StoreHC.Path	d:/dicom/jdicom/print


You can also use the \ to separate the subdirectories.

The stored files which you can find under the StoreHC.Path can be loaded and viewed in any DICOM viewer (e.g. eFilm, DICOMscope, Osiris ..)


Validate

With the JDICOM application Validate you can check DICOM files according the conformity to the DICOM standard.

The JDICOM application Validate runs with a graphical user interface (GUI). Select the entry



in the JDicom menu.



In the JDICOM Validate application you find several tab cards. The tab cards are related to the IOD definitions.

With "**Validate**" only one selected dicom file can be checked.

With "**Validate***" multiple dicom file can be checked.

With "**HTML**" you can export from the validated DICOM file the attributes in an HTML file.

In the log area you will get informations about the validation.

In the window part "**Options**" you can select (mark ✓) optional DICOM modules for your IOD to be checked.

Example for the validation of a MR image. To display the DICOM tags of the selected image the "Attributes" tab card needs to be selected.

CR	Y/N	Tag	Attribute Name	VR	VM	Value
CT	<input checked="" type="checkbox"/>	(0002,0010)	TransferSyntaxUID	UI	1	1.2.840.10008.1.2.1
MR	<input checked="" type="checkbox"/>	(0002,0012)	ImplementationClassUID	UI	1	1.3.12.2.1107.5.2
NM	<input checked="" type="checkbox"/>	(0002,0013)	ImplementationVersionName	SH	1	MREASE_VA11A
PET	<input checked="" type="checkbox"/>	(0008,0005)	SpecificCharacterSet	CS	1	ISO_IR 100
US	<input checked="" type="checkbox"/>	(0008,0008)	ImageType	CS	3	ORIGINAL\PRIMARY\M
US MF	<input checked="" type="checkbox"/>	(0008,0016)	SOPClassUID	UI	1	1.2.840.10008.5.1.4.1.1.4
SC	<input checked="" type="checkbox"/>	(0008,0018)	SOPInstanceUID	UI	1	2.6.200011240845440734000
XA	<input checked="" type="checkbox"/>	(0008,0020)	StudyDate	DA	1	20001122
RF	<input checked="" type="checkbox"/>	(0008,0021)	SeriesDate	DA	1	20001122
St.OV	<input checked="" type="checkbox"/>	(0008,0022)	AcquisitionDate	DA	1	20001122
St.CV	<input checked="" type="checkbox"/>	(0008,0023)	ContentDate	DA	1	20001122
St.LUT	<input checked="" type="checkbox"/>	(0008,0030)	StudyTime	TM	1	075053.000
GSPS	<input checked="" type="checkbox"/>	(0008,0031)	SeriesTime	TM	1	075218.000
Attributes	<input checked="" type="checkbox"/>	(0008,0032)	AcquisitionTime	TM	1	075552.519
	<input checked="" type="checkbox"/>	(0008,0033)	ContentTime	TM	1	074758.000
	<input checked="" type="checkbox"/>	(0008,0050)	AccessionNumber	SH	0	
	<input checked="" type="checkbox"/>	(0008,0060)	Modality	CS	1	MR
	<input checked="" type="checkbox"/>	(0008,0070)	Manufacturer	LO	1	SIEMENS
	<input checked="" type="checkbox"/>	(0008,0090)	ReferringPhysiciansName	PN	0	
	<input checked="" type="checkbox"/>	(0008,1030)	StudyDescription	LO	1	Sprunggelenk^Routine
	<input checked="" type="checkbox"/>	(0008,103E)	SeriesDescription	LO	1	localizer
	<input checked="" type="checkbox"/>	(0008,1090)	ManufacturerModelName	LO	1	Symphony
	<input checked="" type="checkbox"/>	(0008,2111)	DerivationDescription	ST	1	Force Anonymity
	<input checked="" type="checkbox"/>	(0010,0010)	PatientName	PN	1	Anonymous
	<input checked="" type="checkbox"/>	(0010,0020)	PatientID	LO	1	Anonymous
	<input checked="" type="checkbox"/>	(0010,0030)	PatientBirthDate	DA	0	
	<input checked="" type="checkbox"/>	(0010,0040)	PatientSex	CS	1	M
	<input checked="" type="checkbox"/>	(0010,1010)	PatientAge	AS	1	048Y
	<input checked="" type="checkbox"/>	(0010,1030)	PatientWeight	DS	1	99.0

Validate file D:\Gartner\DICOM\mr22076\mr22076_5.0.9133179.img
 No error detected for MR Image Storage SOP Class


In the log you get the following information :


```
Validate file D:\Gartner\DICOM\mr22076\mr22076_5.0.9133179.img
No error detected for MR Image Storage SOP Class
```


This indicates that the DICOM file was successfully checked.

QueryRetrieveSCU

With the JDICOM application QueryRetrieveSCU you can simulate the behaviour of the DICOM service Query Retrieve as Service Class User. An Association can be opened, a DICOM Verification can be made. The image database of a QueryRetrieveSCP can be queried and images can be retrieved.



The DICOM application QueryRetrieveSCU runs with a graphical user interface (GUI). Therefore select the entry  QueryRetrieveSCU in the JDicom menu.



This are the minimum entries which should be modified before you can use QueryRetrieveSCU:

Host :

Enter here the hostname or IP address of the QueryRetrieveSCP (DICOM partner for querying the database and retrieving images)

If a host name is used, an entry in the hosts file is required or DNS has to be enabled.

Port :

Port number of the QueryRetrieveSCP (DICOM partner querying the database and retrieving images)

Called Title:

DICOM AET of the QueryRetrieveSCP (querying the database and retrieving images)

Calling Title:

Own DICOM AET

MoveDestination:

DICOM AET of the StoreSCP which should receive the retrieved images

InfoModel :

Used information model. Possible entries are:

PatientRoot

StudyRoot

PatientStudyOnly

Verbose :

Specifies the log level. The value which can be used starts with

"0" no output in the log tab window

"1"

"2" decode of DIMSE (DICOM Message Service Element)

"3" decode of data elements

Note: If you change an entry in the Props tab you have to confirm it by pressing the <ENTER> key on your keyboard.

If you change a property value you have to CLOSE/OPEN the association in order to activate the made changes.

In this example the DICOM partner for receiving the Query/Retrieve requests (Q/R SCP) has the following parameters:

Hostname: **146.254.123.173**

Port number of the QueryRetrieveSCP process: **104**

AET of the QueryRetrieveSCP process: **AN_TDSIM1**

MoveDestination AET of the StoreSCP process: **QR_SCU**

(AET of the JDICOM ImageServer on the local PC)

The local (own) DICOM AET of the QueryRetrieveSCU JDICOM program: **Storage_SCP**

Used QueryRetrieve information model: **PatientRoot**

Verbose mode: **1**


Start the Association by clicking the "**Connect**" button.

With "**Echo**" you start a Verification (C_Echo).

With "**Query**" sends an query (C_Find) request to the QueryRetrieveSCP.

With "**Retrieve**" you can retrieve the selected Patient/Study/images according the used information model.


- Filter tab card of the JDICOM QueryRetrieveSCU application




In this filter tabcard you can select the tags which shall be used in the query request. Characters entered in under Value will be sent in the Query request as searching (matching) criteria. You can do this for the DICOM tag groups Patient, Study, Series and Image by selecting the tab card for the corresponding tag group.

- Result tab card of the JDICOM QueryRetrieveSCU application


After sending an query request (C_Find) to the QueryRetrieveSCP the feedback is displayed in the result tab card. In this example the information model **PatientRoot** was used.



Selecting the Patient name in the Query Result tree you get DICOM tags displayed for the patient. To query the database for study information you can do this either by double clicking on the patient name or by pressing the "Query" button . The study level information are than visible.



Selecting the study (folder icon) you get then the study DICOM tags displayed. In order to get the information of the image level you can do this either by double clicking on the study folder or by pressing the "Query" button. The image level information are than visible.





Selecting the image in the Query tree you get the image DICOM tags displayed.

You can retrieve the images by selecting in

- ◆ the patient level the patient name and click the "**Retrieve**" button.
- ◆ the study level the study folder and click the "**Retrieve**" button.
- ◆ the image level the image and click the "**Retrieve**" button.



Be sure that the StoreSCP (e.g. JDICOM ImageServer) is running, otherwise the retrieve will fail.

With  you can remove the selected entry in the Query Results tree and
with  you can clear the Query Results tree.

EditDicomObjects

With the JDICOM application EditDicomObject you can create new or modify DICOM files. This files can be DICOM images, worklist entries or mpps files.

The JDICOM application EditDicomObjects runs with a graphical user interface (GUI). Select the entry



Description of the following buttons:



With the **"New"** button you can create a new DICOM file.

With **"Load"** you can load an existing DICOM file (image file or worklist file from JDICOM WL)

With **"Save"** you can save the modified file.

With **"HTML"** you can create a dump of the used DICOM tags from the loaded DICOM file in HTML format.

Note:

Check value of the Property "Save . FileMetaInfo".

If the Value entry of this property is "true" you have to fill out the following File Meta Information Tag:


Tag	Attribute Name	VR	VM	Value
(0002,0012)	ImplementationClassUID	UI	1	nnn.nnn.nnn



If this entry isn't made you will get the error message "com.archimed.dicom.DicomException: Missing ImplementationClassUID" when you try to save the file. In order to avoid this you can modify the value entry in the property "Save.FileMetaInfo" to **false**. This change can not be saved from the GUI EditDicomObjects. When you start the JDicom program EditDicomObjects the value for the tag "Save.FileMetaInfo" is set back to the original one. But there is a way to modify the tag by changing the value in the EditDicomObjects.properties file. You will find this file in the installation directory of JDicom. Use a text editor (e.g. notepad) in order to do the necessary modification. Contents of the "EditDicomObjects.properties" file.


```
#Properties for EditDicomObject
#Thu Jul 19 16:22:40 GMT+02:00 2001
Save.UndefSeqLen=false
Load.PixelData=true
Save.FileMetaInfo=true          (property which should be set to false)
Load.Append=false
Save.TransferSyntax=ImplicitVRLittleEndian
Show.FileMetaInfo=true
Show.SequenceItem=true
Verbose=1
```

Procedure how to insert or remove tags :

- ◆ Load the DICOM file into the EditDicomObjects editor
- ◆ Choose from the IOD's window the part which you want to insert




For inserting the tags select the entry in the IOD's tree with the mouse and use the  button to add the tags. If you want to remove the tags, select first the entry in the IOD's tree and use the  to remove the selected IOD's tag group.



By double clicking on the folder icon of the IOD's tree you can expand the tree for the required selection (e.g. IOD's for MR). Use again the buttons



To add or to remove MR IOD's from the loaded DICOM file.

For example if you want to add the group Patient from the IOD's MR you have to expand the folder MR in the IOD's tree select the entry "Patient" and use the add button  to get the Patient MR IOD tags. You will get than the following output:

Y/N	Tag	Attribute Name	VR	VM	Value
<input type="checkbox"/>	(0008,1120)	ReferencedPatientSequence	SQ	-	
<input checked="" type="checkbox"/>	(0010,0010)	PatientName	PN	1	Hallo^Hallo
<input type="checkbox"/>	(0010,0020)	PatientID	LO	-	
<input type="checkbox"/>	(0010,0030)	PatientBirthDate	DA	-	
<input type="checkbox"/>	(0010,0032)	PatientBirthTime	TM	-	
<input type="checkbox"/>	(0010,0040)	PatientSex	CS	-	
<input type="checkbox"/>	(0010,1000)	OtherPatientID	LO	-	
<input type="checkbox"/>	(0010,1001)	OtherPatientNames	PN	-	
<input type="checkbox"/>	(0010,2160)	EthnicGroup	SH	-	
<input type="checkbox"/>	(0010,4000)	PatientComments	LT	-	

Select the tag which you want to use and insert the correct value. Automatically the tag gets marked with ✓ after confirming the made entry with the <ENTER> key on the keyboard. For resetting the value remove the mark ✓ by clicking on it.


Note: If you change or make an IOD entry you have to confirm it by pressing the <ENTER> key on your keyboard.

If all changes were made on the DICOM file, save it.

EditDicomDir

The JDICOM application EditDicomDir runs with a graphical user interface (GUI). Select the entry

 EditDicomDir in the JDicom menu.



Description of the buttons:

New

with the new button an empty DICOMDIR is generated, but not save to disc.

Load

with the load button you can load an existing DICOMDIR file into the JDICOM application EditDicomDir.

Save

with the save button you can save the DICOMDIR file to disc.

Example for a loaded DICOMDIR file:

